

My Peristeen[®] Plus irrigation checklist

Good preparation is key when performing transanal irrigation. Make sure you have everything you are going to need before you begin. Use this checklist to set yourself up for success.

(Print me out or stick me to your bathroom door for an easy reminder before every irrigation!)

You need:

- Your Peristeen Plus system + Instructions for Use
- Wipes
- Disposal bag

You may want:

- Leg straps
- Disposable gloves
- Paper towels
- Something to entertain yourself during your waiting time
- A phone, in case you wish to contact someone for assistance

Now you have everything you need, have a quick check – is there anything that may prevent you from completing your irrigation as you would like?

Always make sure you have enough time to complete your irrigation. Everyone is different, but it takes on average 30 minutes to empty the bowels after pumping the water.

I have sufficient time to complete my irrigation

Peristeen Plus users often wish to complete their irrigations in a private and quiet setting, without being disturbed. If you live with others, you may wish to let them know you need some privacy.

I am in a comfortable environment where I won't be disturbed

You are ready to get started.

Don't forget you can fill out your irrigation diary after you complete your routine to track your progress.



Scan to get your
irrigation diary or visit
coloplastcare.ca/bowel

And remember, your Coloplast Care Advisor is available (1-866-293-6349, 9:00AM -5:00PM ET) if you require further support.

[Ostomy Care](#) | [Continence Care](#) | [Wound and Skin Care](#) | [Interventional Urology](#) | [Voice and Respiratory Care](#)

Information from Coloplast Care is for educational purposes only. It is not intended to substitute for professional medical advice and should not be interpreted to contain treatment recommendations. Talk to your healthcare professional about what may be right for you.

Coloplast Canada: 205A - 2401 Bristol Circle, Oakville, ON, Canada 1-866-293-6349

www.coloplast.ca The Coloplast logo is a registered trademark of Coloplast A/S. © 2023-03. All rights reserved Coloplast A/S. PM-24276

